

For Research Use Only

# SHOC2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84917-4-PBS



## Basic Information

Catalog Number:

84917-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11754

GenBank Accession Number:

BC050445

GeneID (NCBI):

8036

UNIPROT ID:

Q9UQ13

Full Name:

soc-2 suppressor of clear homolog (C. elegans)

Calculated MW:

582 aa, 65 kDa

Purification Method:

Protein A purification

CloneNo.:

242343F2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01674-2, SHOC2/Sur-8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84917-2-PBS. Detection antibody: 84917-4-PBS. Standard: Ag11754. Range: 0.781-100 ng/mL