

For Research Use Only

PSMD3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84928-4-PBS



Basic Information

Catalog Number: 84928-4-PBS	GenBank Accession Number: BC012302	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5709	CloneNo.: 242470F5
Source: Rabbit	UNIPROT ID: O43242	
Isotype: IgG	Full Name: proteasome (prosome, macropain) 26S subunit, non-ATPase, 3	
Immunogen Catalog Number: AG2684	Calculated MW: 534 aa, 61 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01691-1, PSMD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84928-3-PBS. Detection antibody: 84928-4-PBS. Standard: Ag2684. Range: 0.781-100 ng/mL.