

For Research Use Only

# TYRO3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84929-4-PBS



## Basic Information

Catalog Number: 84929-4-PBS	GenBank Accession Number: BC049368	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7301	CloneNo.: 242477E1
Source: Rabbit	UNIPROT ID: Q06418	
Isotype: IgG	Full Name: TYRO3 protein tyrosine kinase	
Immunogen Catalog Number: AG29527	Calculated MW: 97 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP01701-3, Human TYRO3 Recombinant Matched Antibody Pair - PBS only. 84929-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag29527. 84929-5-PBS was HRP conjugated as the detection antibody. Range: 0.781-50 ng/mL.



The mean TYRO3 concentration was determined to be 19.82 ng/mL in HT-29 cell extract based on a 1.50 mg/mL extract load, 25.19 ng/mL in COLO 320 cell extract based on a 1.40 mg/mL extract load and 20.61 ng/mL in SGC-7901 cell extract based on a 1.60 mg/mL extract load.