

For Research Use Only

TYRO3 Recombinant antibody, PBS Only

Catalog Number: 84929-6-PBS



Basic Information

Catalog Number: 84929-6-PBS	GenBank Accession Number: BC049368	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7301	CloneNo.: 242477G1
Source: Rabbit	UNIPROT ID: Q06418	
Isotype: IgG	Full Name: TYRO3 protein tyrosine kinase	
Immunogen Catalog Number: AG29527	Calculated MW: 97 kDa	
	Observed MW: 130 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84929-6-RR (TYRO3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84929-6-PBS in a different storage buffer formulation.