

For Research Use Only

HGD Recombinant antibody, PBS Only (Capture)

Catalog Number: 84931-3-PBS



Basic Information

Catalog Number: 84931-3-PBS	GenBank Accession Number: BC020792	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3081	CloneNo.: 242485D8
Source: Rabbit	UNIPROT ID: Q93099	
Isotype: IgG	Full Name: homogentisate 1,2-dioxygenase (homogentisate oxidase)	
Immunogen Catalog Number: AG9544	Calculated MW: 37 kDa, 50 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01681-2, HGD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84931-3-PBS. Detection antibody: 84931-1-PBS. Standard: Ag9544. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01681-1, HGD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84931-3-PBS. Detection antibody: 84931-2-PBS. Standard: Ag9544. Range: 0.625-80 ng/mL