

For Research Use Only

YOD1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84933-3-PBS



Basic Information

Catalog Number: 84933-3-PBS	GenBank Accession Number: BC137167	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55432	CloneNo.: 242498D5
Source: Rabbit	UNIPROT ID: Q5VVQ6	
Isotype: IgG	Full Name: YOD1 OTU deubiquinating enzyme 1 homolog (S. cerevisiae)	
Immunogen Catalog Number: AG21844	Calculated MW: 348 aa, 38 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01684-1, YOD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84933-4-PBS. Detection antibody: 84933-3-PBS. Standard: Ag21844. Range: 0.781-100 ng/mL.