

For Research Use Only

# IL-5 R alpha/CD125 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84948-3-PBS



## Basic Information

Catalog Number:

84948-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_175726.4

GeneID (NCBI):

3568

UNIPROT ID:

Q01344-1

Full Name:

interleukin 5 receptor, alpha

Calculated MW:

48kDa

Purification Method:

Protein A purification

CloneNo.:

242234H6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01699-2, IL-5RA/CD125 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84948-3-PBS. Detection antibody: 84948-2-PBS. Standard: Eg1665. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01699-1, IL-5RA/CD125 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84948-1-PBS. Detection antibody: 84948-3-PBS. Standard: Eg1665. Range: 0.625-80 ng/mL