

For Research Use Only

HERC4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84952-3-PBS



Basic Information

Catalog Number: 84952-3-PBS	GenBank Accession Number: BC039600	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 26091	CloneNo.: 242599G10
Source: Rabbit	UNIPROT ID: Q5GLZ8	
Isotype: IgG	Full Name: hct domain and RLD 4	
Immunogen Catalog Number: AG4566	Calculated MW: 1049 aa, 118 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01706-1, HERC4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84952-3-PBS. Detection antibody: 84952-4-PBS. Standard: Ag4566. Range: 0.313-40 ng/mL