

For Research Use Only

CSTA Recombinant antibody, PBS Only (Capture)

Catalog Number: 84953-3-PBS



Basic Information

Catalog Number: 84953-3-PBS	GenBank Accession Number: BC010379	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1475	CloneNo.: 242558G8
Source: Rabbit	UNIPROT ID: P01040	
Isotype: IgG	Full Name: cystatin A (stefin A)	
Immunogen Catalog Number: AG8733	Calculated MW: 98 aa, 11 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01707-2, Cystatin A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84953-3-PBS. Detection antibody: 84953-1-PBS. Standard: Ag8733. Range: 0.469-60 ng/mL.