

For Research Use Only

TFB1M Recombinant antibody, PBS Only

Catalog Number: 84957-4-PBS



Basic Information

Catalog Number: 84957-4-PBS	GenBank Accession Number: BC017788	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 51106	CloneNo.: 242524A10
Source: Rabbit	UNIPROT ID: Q8WVMO	
Isotype: IgG	Full Name: transcription factor B1, mitochondrial	
Immunogen Catalog Number: AG9785	Calculated MW: 346 aa, 40 kDa	
	Observed MW: 40 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84957-4-RR (TFB1M antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84957-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84957-4-RR against Human TFB1M were performed. The affinity constant is 2.99 nM.