

For Research Use Only

SUGT1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84963-1-PBS



Basic Information

Catalog Number: 84963-1-PBS	GenBank Accession Number: BC000911	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10910	CloneNo.: 242590C9
Source: Rabbit	UNIPROT ID: Q9Y2ZO	
Isotype: IgG	Full Name: SGT1, suppressor of G2 allele of SKP1 (S. cerevisiae)	
Immunogen Catalog Number: AG2275	Calculated MW: 41 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01689-2, SUGT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84963-2-PBS. Detection antibody: 84963-1-PBS. Standard: Ag2275. Range: 0.469-60 ng/mL.