

For Research Use Only

# Eg5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84966-1-PBS



## Basic Information

Catalog Number: 84966-1-PBS	GenBank Accession Number: BC136474	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3832	CloneNo.: 242611C10
Source: Rabbit	UNIPROT ID: P52732	
Isotype: IgG	Full Name: kinesin family member 11	
Immunogen Catalog Number: AG20001	Calculated MW: 1056 aa, 119 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01688-2, Eg5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84966-2-PBS. Detection antibody: 84966-1-PBS. Standard: Ag20001. Range: 3.125-100 ng/mL.