

For Research Use Only

MAP2K2 Recombinant antibody, PBS Only (Detector)

Catalog Number:84967-3-PBS



Basic Information

Catalog Number: 84967-3-PBS	GenBank Accession Number: BC018645	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5605	CloneNo.: 242614H4
Source: Rabbit	UNIPROT ID: P36507	
Isotype: IgG	Full Name: mitogen-activated protein kinase kinase 2	
Immunogen Catalog Number: AG1481	Calculated MW: 44 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01692-2, MEK1/2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84967-1-PBS. Detection antibody: 84967-3-PBS. Standard: Ag1481. Range: 0.625-80 ng/mL