

For Research Use Only

GFER Recombinant antibody, PBS Only (Capture)

Catalog Number: 84972-1-PBS



Basic Information

Catalog Number: 84972-1-PBS	GenBank Accession Number: BC028348	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2671	CloneNo.: 242557E2
Source: Rabbit	UNIPROT ID: P55789	
Isotype: IgG	Full Name: growth factor, augmenter of liver regeneration	
Immunogen Catalog Number: AG1840	Calculated MW: 15 kDa, 23 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01721-1, ALR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84972-1-PBS. Detection antibody: 84972-3-PBS. Standard: Ag1840. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP01721-2, ALR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84972-1-PBS. Detection antibody: 84972-2-PBS. Standard: Ag1840. Range: 0.156-20 ng/mL