

For Research Use Only

C3orf20 Recombinant antibody, PBS Only

Catalog Number: 84975-1-PBS



Basic Information

Catalog Number: 84975-1-PBS	GenBank Accession Number: BC038406	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 84077	CloneNo.: 242404E2
Source: Rabbit	UNIPROT ID: Q8ND61	
Isotype: IgG	Full Name: chromosome 3 open reading frame 20	
Immunogen Catalog Number: AG15031	Calculated MW: 904 aa, 101 kDa	
	Observed MW: 88 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

The function of C3orf20 has not been widely studied, and is yet to be fully elucidated. Two isoforms were produced by alternative splicing with predicted MW of 88 and 101 kDa according to UniProt. Catalog#20906-1-AP can recognize both 88 and 101 kDa isoforms.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84975-1-RR (C3orf20 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84975-1-PBS in a different storage buffer formulation.