

For Research Use Only

TRPA1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84986-1-PBS



Basic Information

Catalog Number: 84986-1-PBS	GenBank Accession Number: NM_007332	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8989	CloneNo.: 242162A2
Source: Rabbit	UNIPROT ID: O75762	
Isotype: IgG	Full Name: transient receptor potential cation channel, subfamily A, member 1	
	Calculated MW: 140 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01728-1, Human TRPA1 Recombinant Matched Antibody Pair - PBS only. 84986-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY00327. 84986-2-PBS was HRP conjugated as the detection antibody. Range: 0.312-20 ng/mL

The mean TRPA1 concentration was determined to be 331.00 ng/mL in SH-SY5Y cell extract based on a 1.0 mg/mL extract load.