

For Research Use Only

TRPA1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84986-2-PBS



Basic Information

Catalog Number:

84986-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_007332

GeneID (NCBI):

8989

UNIPROT ID:

O75762

Full Name:

transient receptor potential cation
channel, subfamily A, member 1

Calculated MW:

140 kDa

Purification Method:

Protein A purification

CloneNo.:

242162D11

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01728-1, Human TRPA1 Recombinant Matched Antibody Pair - PBS only. 84986-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY00327. 84986-2-PBS was HRP conjugated as the detection antibody. Range: 0.312-20 ng/mL

The mean TRPA1 concentration was determined to be 331.00 ng/mL in SH-SY5Y cell extract based on a 1.0 mg/mL extract load.