

For Research Use Only

IL-3 R Beta/CD131 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84988-2-PBS



Basic Information

Catalog Number: 84988-2-PBS	GenBank Accession Number: NM_000395.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1439	CloneNo.: 241823B6
Source: Rabbit	UNIPROT ID: P32927-1	
Isotype: IgG	Full Name: colony stimulating factor 2 receptor, beta, low-affinity (granulocyte-macrophage)	
	Calculated MW: 97kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01730-1, IL-3 R Beta/CD131 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84988-2-PBS. Detection antibody: 84988-4-PBS. Standard: Eg1662. Range: 0.625-20 ng/mL.