

For Research Use Only

# TRMT10C Recombinant antibody, PBS Only (Detector)

Catalog Number: 84991-2-PBS



## Basic Information

Catalog Number:

84991-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30642

GenBank Accession Number:

NM\_017819

GeneID (NCBI):

54931

UNIPROT ID:

Q7L0Y3

Full Name:

RNA (guanine-9-) methyltransferase domain containing 1

Calculated MW:

47 kDa

Purification Method:

Protein A purification

CloneNo.:

242522B7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01765-2, TRMT10C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84991-3-PBS. Detection antibody: 84991-2-PBS. Standard: Ag30642. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP01765-1, TRMT10C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84991-1-PBS. Detection antibody: 84991-2-PBS. Standard: Ag30642. Range: 0.625-80 ng/mL.