

For Research Use Only

SNRNP200 Recombinant antibody, PBS Only (Capture)

Catalog Number:84992-1-PBS



Basic Information

Catalog Number: 84992-1-PBS	GenBank Accession Number: BC001417	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 23020	CloneNo.: 242550A7
Source: Rabbit	UNIPROT ID: O75643	
Isotype: IgG	Full Name: small nuclear ribonucleoprotein	
Immunogen Catalog Number: AG20748	200kDa (U5)	
	Calculated MW: 494 aa, 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01734-1, SNRNP200 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84992-1-PBS. Detection antibody: 84992-3-PBS. Standard: Ag20748. Range: 1.563-100 ng/mL.