

For Research Use Only

SLC27A6 Recombinant antibody, PBS Only (Detector)

Catalog Number:84995-3-PBS



Basic Information

Catalog Number:

84995-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG32179

GenBank Accession Number:

BC041945

GeneID (NCBI):

28965

UNIPROT ID:

Q9Y2P4

Full Name:

solute carrier family 27 (fatty acid transporter), member 6

Calculated MW:

70 kDa

Purification Method:

Protein A purification

CloneNo.:

242542H4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01755-1, SLC27A6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84995-1-PBS. Detection antibody: 84995-3-PBS. Standard: Ag32179. Range: 0.781-100 ng/mL.