

For Research Use Only

# NADK Recombinant antibody, PBS Only (Capture)

Catalog Number: 84996-3-PBS



## Basic Information

Catalog Number: 84996-3-PBS	GenBank Accession Number: BC001709	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 65220	CloneNo.: 242502G3
Source: Rabbit	UNIPROT ID: O95544	
Isotype: IgG	Full Name: NAD kinase	
Immunogen Catalog Number: AG7906	Calculated MW: 49 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP01732-2, Human NADK Recombinant Matched Antibody Pair - PBS only. 84996-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag7906. 84996-2-PBS was HRP conjugated as the detection antibody. Range: 0.625-40 ng/mL

The mean NADK concentration was determined to be 26.13 ng/mL in HeLa cell extract based on a 1.5 mg/mL extract load.