

For Research Use Only

TBX5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84997-3-PBS



Basic Information

Catalog Number: 84997-3-PBS	GenBank Accession Number: BC027942	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6910	CloneNo.: 242565D9
Source: Rabbit	UNIPROT ID: Q99593	
Isotype: IgG	Full Name: T-box 5	
Immunogen Catalog Number: AG3845	Calculated MW: 518 aa, 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01748-1, TBX5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84997-2-PBS. Detection antibody: 84997-3-PBS. Standard: Ag3845. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01748-2, TBX5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84997-1-PBS. Detection antibody: 84997-3-PBS. Standard: Ag3845. Range: 1.563-100 ng/mL.