

For Research Use Only

ZNF618 Recombinant antibody, PBS Only (Capture)

Catalog Number:84998-3-PBS



Basic Information

Catalog Number: 84998-3-PBS	GenBank Accession Number: NM_001318042.1	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 114991	CloneNo.: 242566G9
Source: Rabbit	UNIPROT ID: Q5T7W0	
Isotype: IgG	Full Name: zinc finger protein 618	
Immunogen Catalog Number: AG36487	Calculated MW: 105kDa, 954aa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01753-1, ZNF618 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84998-3-PBS. Detection antibody: 84998-4-PBS. Standard: Ag36487. Range: 0.313-40 ng/mL.