

For Research Use Only

DCBLD2 Recombinant antibody, PBS Only (Capture)

Catalog Number:85000-1-PBS



Basic Information

Catalog Number: 85000-1-PBS	GenBank Accession Number: BC029658	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 131566	CloneNo.: 242604B6
Source: Rabbit	UNIPROT ID: Q96PD2	
Isotype: IgG	Full Name: discoidin, CUB and LCCL domain containing 2	
Immunogen Catalog Number: AG32823	Calculated MW: 775 aa, 89 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01733-2, DCBLD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85000-1-PBS. Detection antibody: 85000-2-PBS. Standard: Ag32823. Range: 0.469-15 ng/mL.