

For Research Use Only

# MAS1L Recombinant antibody

Catalog Number:85003-1-RR



## Basic Information

Catalog Number:

85003-1-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_052967

GeneID (NCBI):

116511

UNIPROT ID:

P35410

Full Name:

MAS1 oncogene-like

Calculated MW:

42 kDa

Observed MW:

42 kDa

Purification Method:

Protein A purification

CloneNo.:

242360D1

Recommended Dilutions:

WB 1:2000-1:10000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HeLa cells, HepG2 cells, HEK-293 cells, COLO 320 cells, NIH/3T3 cells

## Background Information

MAS1L, also named MRG, belongs to the G-protein coupled receptor 1 family and Mas subfamily.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85003-1-RR (MAS1L antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.