

For Research Use Only

LRRC41 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85012-1-PBS



Basic Information

Catalog Number:

85012-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14266

GenBank Accession Number:

BC004953

GeneID (NCBI):

10489

UNIPROT ID:

Q15345

Full Name:

leucine rich repeat containing 41

Calculated MW:

812 aa, 87 kDa

Purification Method:

Protein A purification

CloneNo.:

242380F10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01742-1, LRRC41 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85012-2-PBS. Detection antibody: 85012-1-PBS. Standard: Ag14266. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01742-2, LRRC41 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85012-2-PBS. Detection antibody: 85012-3-PBS. Standard: Ag14266. Range: 0.781-100 ng/mL.