

For Research Use Only

# LRRC41 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85012-2-PBS



## Basic Information

Catalog Number: 85012-2-PBS	GenBank Accession Number: BC004953	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10489	CloneNo.: 242380F2
Source: Rabbit	UNIPROT ID: Q15345	
Isotype: IgG	Full Name: leucine rich repeat containing 41	
Immunogen Catalog Number: AG14266	Calculated MW: 812 aa, 87 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01742-1, LRRC41 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85012-2-PBS. Detection antibody: 85012-1-PBS. Standard: Ag14266. Range: 0.781-100 ng/mL.