

For Research Use Only

CPT2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85013-3-PBS



Basic Information

Catalog Number: 85013-3-PBS	GenBank Accession Number: BC002445	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1376	CloneNo.: 242381G9
Source: Rabbit	UNIPROT ID: P23786	
Isotype: IgG	Full Name: carnitine palmitoyltransferase 2	
Immunogen Catalog Number: AG24897	Calculated MW: 74 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01757-1, CPT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85013-3-PBS. Detection antibody: 85013-2-PBS. Standard: Ag24897. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01757-2, CPT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85013-3-PBS. Detection antibody: 85013-1-PBS. Standard: Ag24897. Range: 0.781-100 ng/mL.