

For Research Use Only

Selenoprotein P Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number:85017-2-PBS



Basic Information

Catalog Number: 85017-2-PBS	GenBank Accession Number: NM_001085486.2	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6414	CloneNo.: 242510D6
Source: Rabbit	UNIPROT ID: P49908	
Isotype: IgG	Full Name: selenoprotein P, plasma, 1	
Immunogen Catalog Number: AG32952	Calculated MW: 43 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01743-2, Selenoprotein P Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85017-3-PBS. Detection antibody: 85017-2-PBS. Standard: Ag32952. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01743-1, Selenoprotein P Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85017-2-PBS. Detection antibody: 85017-1-PBS. Standard: Ag32952. Range: 0.781-100 ng/mL