

For Research Use Only

SMARCA5 Recombinant antibody, PBS Only (Detector)

Catalog Number:85019-1-PBS



Basic Information

Catalog Number:

85019-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3759

GenBank Accession Number:

BC023144

GeneID (NCBI):

8467

UNIPROT ID:

O60264

Full Name:

SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 5

Calculated MW:

1052 aa, 122 kDa

Purification Method:

Protein A purification

CloneNo.:

242512A10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01767-2, SMARCA5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85019-3-PBS. Detection antibody: 85019-1-PBS. Standard: Ag3759. Range: 0.781-100 ng/mL.