

For Research Use Only

# SMARCA5 Recombinant antibody, PBS Only (Detector)

Catalog Number:85019-4-PBS



## Basic Information

Catalog Number:

85019-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3759

GenBank Accession Number:

BC023144

GeneID (NCBI):

8467

UNIPROT ID:

O60264

Full Name:

SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 5

Calculated MW:

1052 aa, 122 kDa

Purification Method:

Protein A purification

CloneNo.:

242512D11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP0176-1, SMARCA5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85019-2-PBS. Detection antibody: 85019-4-PBS. Standard: Ag3759. Range: 0.781-100 ng/mL.