

For Research Use Only

# RSPH9 Recombinant antibody, PBS Only (Capture)

Catalog Number:85026-1-PBS



## Basic Information

<b>Catalog Number:</b> 85026-1-PBS	<b>GenBank Accession Number:</b> BC029519	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 221421	<b>CloneNo.:</b> 24258483
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9H1X1	
<b>Isotype:</b> IgG	<b>Full Name:</b> radial spoke head 9 homolog (Chlamydomonas)	
<b>Immunogen Catalog Number:</b> AG19745	<b>Calculated MW:</b> 276 aa, 31 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01764-2, RSPH9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85026-1-PBS. Detection antibody: 85026-3-PBS. Standard: Ag19745. Range: 0.625-40 ng/mL