

For Research Use Only

PDLIM1,CLP36 Recombinant antibody, PBS Only (Detector)

Catalog Number:85027-1-PBS



Basic Information

Catalog Number: 85027-1-PBS	GenBank Accession Number: BC000915	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9124	CloneNo.: 242585C8
Source: Rabbit	UNIPROT ID: O00151	
Isotype: IgG	Full Name: PDZ and LIM domain 1	
Immunogen Catalog Number: AG2274	Calculated MW: 329 aa, 36 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01749-1, PDLIM1/CLP36 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85027-3-PBS. Detection antibody: 85027-1-PBS. Standard: Ag2274. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01749-2, PDLIM1/CLP36 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85027-2-PBS. Detection antibody: 85027-1-PBS. Standard: Ag2274. Range: 0.313-40 ng/mL