

For Research Use Only

# PDLIM1,CLP36 Recombinant antibody, PBS Only (Capture)

Catalog Number:85027-2-PBS



## Basic Information

Catalog Number: 85027-2-PBS	GenBank Accession Number: BC000915	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9124	CloneNo.: 242585E4
Source: Rabbit	UNIPROT ID: O00151	
Isotype: IgG	Full Name: PDZ and LIM domain 1	
Immunogen Catalog Number: AG2274	Calculated MW: 329 aa, 36 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01749-2, PDLIM1/CLP36 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85027-2-PBS. Detection antibody: 85027-1-PBS. Standard: Ag2274. Range: 0.313-40 ng/mL.