

For Research Use Only

HEXIM1 Recombinant antibody, PBS Only (Capture)

Catalog Number:85028-3-PBS



Basic Information

Catalog Number: 85028-3-PBS	GenBank Accession Number: BC006460	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10614	CloneNo.: 242587F1
Source: Rabbit	UNIPROT ID: O94992	
Isotype: IgG	Full Name: hexamethylene bis-acetamide inducible 1	
Immunogen Catalog Number: AG8144	Calculated MW: 41 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01740-2, HEXIM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85028-3-PBS. Detection antibody: 85028-1-PBS. Standard: Ag8144. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP01740-1, HEXIM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85028-3-PBS. Detection antibody: 85028-2-PBS. Standard: Ag8144. Range: 0.313-20 ng/mL