

For Research Use Only

C1R Recombinant antibody, PBS Only (Capture)

Catalog Number:85029-2-PBS



Basic Information

Catalog Number: 85029-2-PBS	GenBank Accession Number: BC035220	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 715	CloneNo.: 242596E2
Source: Rabbit	UNIPROT ID: P00736	
Isotype: IgG	Full Name: complement component 1, r subcomponent	
Immunogen Catalog Number: AG11390	Calculated MW: 705 aa, 80 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01761-2, C1r Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85029-2-PBS. Detection antibody: 85029-3-PBS. Standard: Ag11390. Range: 0.781-100 ng/mL.