

For Research Use Only

# C1QC Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number:85031-1-PBS



## Basic Information

Catalog Number:

85031-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10400

GenBank Accession Number:

BC009016

GeneID (NCBI):

714

UNIPROT ID:

P02747

Full Name:

complement component 1, q  
subcomponent, C chain

Calculated MW:

245 aa, 26 kDa

Purification Method:

Protein A purification

CloneNo.:

242624A2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01746-1, C1qC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85031-1-PBS. Detection antibody: 85031-2-PBS. Standard: Ag10400. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01746-2, C1qC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85031-3-PBS. Detection antibody: 85031-1-PBS. Standard: Ag10400. Range: 0.781-100 ng/mL.