

For Research Use Only

SEC22B Recombinant antibody, PBS Only (Detector)

Catalog Number:85037-2-PBS



Basic Information

Catalog Number: 85037-2-PBS	GenBank Accession Number: BC001364	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9554	CloneNo.: 242659C7
Source: Rabbit	UNIPROT ID: O75396	
Isotype: IgG	Full Name: SEC22 vesicle trafficking protein homolog B (S. cerevisiae)	
Immunogen Catalog Number: AG12981	Calculated MW: 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01774-2, SEC22B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85037-1-PBS. Detection antibody: 85037-2-PBS. Standard: Ag12981. Range: 0.313-40 ng/mL