

For Research Use Only

DCP1A Recombinant antibody, PBS Only (Detector)

Catalog Number: 85047-3-PBS



Basic Information

Catalog Number: 85047-3-PBS	GenBank Accession Number: BC007439	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55802	CloneNo.: 242553E12
Source: Rabbit	UNIPROT ID: Q9NPI6	
Isotype: IgG	Full Name: DCP1 decapping enzyme homolog A (S. cerevisiae)	
Immunogen Catalog Number: AG18027	Calculated MW: 582 aa, 63 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01815-1, DCP1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85047-2-PBS. Detection antibody: 85047-3-PBS. Standard: Ag18027. Range: 0.313-10 ng/mL



Cytometric bead array standard curve of MP01815-2, DCP1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85047-1-PBS. Detection antibody: 85047-3-PBS. Standard: Ag18027. Range: 0.313-10 ng/mL