

For Research Use Only

PML Recombinant antibody, PBS Only (Capture)

Catalog Number:85048-1-PBS



Basic Information

Catalog Number: 85048-1-PBS	GenBank Accession Number: BC020994	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5371	CloneNo.: 242564C12
Source: Rabbit	UNIPROT ID: P29590	
Isotype: IgG	Full Name: promyelocytic leukemia	
Immunogen Catalog Number: AG15377	Calculated MW: 882 aa, 98 kDa	

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA, Indirect ELISA,
Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP01806-2, Human PML Recombinant Matched Antibody Pair - PBS only. 85048-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag15377. 85048-3-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



The mean PML concentration was determined to be 5.94 ng/mL in HeLa cell extract based on a 1.50 mg/mL extract load.



Cytometric bead array standard curve of MP01806-2, PML Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85048-1-PBS. Detection antibody: 85048-3-PBS. Standard: Ag15377. Range: 0.781-100 ng/mL