

For Research Use Only

SP3 Recombinant antibody, PBS Only (Detector)

Catalog Number:85050-3-PBS



Basic Information

Catalog Number: 85050-3-PBS	GenBank Accession Number: BC126414	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6670	CloneNo.: 242556G9
Source: Rabbit	UNIPROT ID: Q02447	
Isotype: IgG	Full Name: Sp3 transcription factor	
Immunogen Catalog Number: AG25082	Calculated MW: 781 aa, 82 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01797-2, SP3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85050-2-PBS. Detection antibody: 85050-3-PBS. Standard: Ag25082. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP01797-1, SP3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85050-1-PBS. Detection antibody: 85050-3-PBS. Standard: Ag25082. Range: 0.156-20 ng/mL