

For Research Use Only

SP3 Recombinant antibody, PBS Only (Detector)

Catalog Number:85050-4-PBS



Basic Information

Catalog Number: 85050-4-PBS	GenBank Accession Number: BC126414	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6670	CloneNo.: 242556A3
Source: Rabbit	UNIPROT ID: Q02447	
Isotype: IgG	Full Name: Sp3 transcription factor	
Immunogen Catalog Number: AG25082	Calculated MW: 781 aa, 82 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01797-3, Human SP3 Recombinant Matched Antibody Pair - PBS only. 85050-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag25082. 85050-4-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL.



The mean SP3 concentration was determined to be 27.89 ng/mL in PC-3 cell extract based on a 2.50 mg/mL extract load and 32.69 ng/mL in Jurkat cell extract based on a 2.00 mg/mL extract load.