

For Research Use Only

Cannabinoid receptor 2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85051-2-PBS



Basic Information

Catalog Number: 85051-2-PBS	GenBank Accession Number: BC074767	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1269	CloneNo.: 242594F12
Source: Rabbit	UNIPROT ID: P34972	
Isotype: IgG	Full Name: cannabinoid receptor 2 (macrophage)	
Immunogen Catalog Number: AG30609	Calculated MW: 360 aa, 40 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01799-2, Cannabinoid receptor 2/CNR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85051-2-PBS. Detection antibody: 85051-3-PBS. Standard: Ag30609. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP01799-1, Cannabinoid receptor 2/CNR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85051-2-PBS. Detection antibody: 85051-1-PBS. Standard: Ag30609. Range: 0.625-80 ng/mL.