

For Research Use Only

RGS12 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85057-1-PBS



Basic Information

Catalog Number: 85057-1-PBS	GenBank Accession Number: NM_198229	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6002	CloneNo.: 242530F1
Source: Rabbit	UNIPROT ID: O14924	
Isotype: IgG	Full Name: regulator of G-protein signaling 12	
Immunogen Catalog Number: AG30731	Calculated MW: 156 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01790-1, RGS12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85057-1-PBS. Detection antibody: 85057-2-PBS. Standard: Ag30731. Range: 0.625-40 ng/mL