

For Research Use Only

# S1PR1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85059-3-PBS



## Basic Information

Catalog Number:

85059-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28083

GenBank Accession Number:

BC018650

GeneID (NCBI):

1901

UNIPROT ID:

P21453

Full Name:

sphingosine-1-phosphate receptor 1

Calculated MW:

43 kDa

Purification Method:

Protein A purification

CloneNo.:

242592E12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01781-2, S1PR1/EDG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85059-3-PBS. Detection antibody: 85059-1-PBS. Standard: Ag28083. Range: 0.313-40 ng/mL.