

For Research Use Only

DDX46 Recombinant antibody, PBS Only (Detector)

Catalog Number:85065-2-PBS



Basic Information

Catalog Number: 85065-2-PBS	GenBank Accession Number: BC012304	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9879	CloneNo.: 242715E1
Source: Rabbit	UNIPROT ID: Q7L014	
Isotype: IgG	Full Name: DEAD (Asp-Glu-Ala-Asp) box polypeptide 46	
Immunogen Catalog Number: AG10342	Calculated MW: 1031 aa, 117 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01807-2, DDX46 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85065-3-PBS. Detection antibody: 85065-2-PBS. Standard: Ag10342. Range: 0.781-100 ng/mL.