

For Research Use Only

ARR3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85067-3-PBS



Basic Information

Catalog Number:

85067-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1580

GenBank Accession Number:

BC012096

GeneID (NCBI):

407

UNIPROT ID:

P36575

Full Name:

arrestin 3, retinal (X-arrestin)

Calculated MW:

43 kDa

Purification Method:

Protein A purification

CloneNo.:

242751E8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01810-2, Arrestin C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85067-3-PBS. Detection antibody: 85067-1-PBS. Standard:Ag1580. Range: 0.781-100 ng/mL