

For Research Use Only

# NME2 Recombinant antibody, PBS Only (Capture)

Catalog Number:85074-3-PBS



## Basic Information

Catalog Number:

85074-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13913

GenBank Accession Number:

BC002476

GeneID (NCBI):

4831

UNIPROT ID:

P22392

Full Name:

non-metastatic cells 2, protein  
(NM23B) expressed in

Calculated MW:

152 aa, 17 kDa

Purification Method:

Protein A purification

CloneNo.:

242698F4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01804-1, NME2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85074-3-PBS. Detection antibody: 85074-2-PBS. Standard: Ag13913. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP01804-2, NME2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85074-3-PBS. Detection antibody: 85074-1-PBS. Standard: Ag13913. Range: 1.563-200 ng/mL.