

For Research Use Only

FLT4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85076-4-PBS



Basic Information

Catalog Number:

85076-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_182925

GeneID (NCBI):

2324

UNIPROT ID:

P35916-1

Full Name:

fms-related tyrosine kinase 4

Calculated MW:

153kd

Purification Method:

Protein A purification

CloneNo.:

242430G9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00009-5, FLT4/VEGFR3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85076-2-PBS. Detection antibody: 85076-4-PBS. Standard: Eg2494. Range: 0.781-100 ng/mL.